



**Computer and Augmentative Communication Mountings**



36014



Adjustable tilt with quick release lever



Quick release lock for easy computer/communication device removal.



With rear lock removed, unit slides out from rear of mount - no Velcro or cables to remove



Close up of #36026 frame mount



Adjustable tabs and mounting position

Close up of #36014 frame mount



**Computer and communication device mount**

Fully adjustable mounting system for virtually any laptop computer or communication device. Front, side, and rear holding points are easily adjusted for device width, depth, and thickness. Quick device removal with adjustable locking bar in rear. Mounts to wheelchair or hospital bed frame with adjustable angle quick release receiver. Custom sizing available for holding points, mounts, etc.

ITEM #	DESCRIPTION	MOUNT STYLE AND SIZE
36010	COMPUTER AND COMMUNICATION DEVICE MOUNT	DUAL ROUND COLLARS - 3/4"
36012	COMPUTER AND COMMUNICATION DEVICE MOUNT	DUAL ROUND COLLARS - 7/8"
36014	COMPUTER AND COMMUNICATION DEVICE MOUNT	DUAL ROUND COLLARS - 1"
36020	COMPUTER AND COMMUNICATION DEVICE MOUNT	RECTANGULAR TUBE - 3/4" X 1-1/2"
36022	COMPUTER AND COMMUNICATION DEVICE MOUNT	RECTANGULAR TUBE - 1" X 2"
36026	COMPUTER AND COMMUNICATION DEVICE MOUNT	FLAT PLATE MOUNT - 3" CENTER TO CENTER OF OVAL HOLES




**Augmentative Communication Mounting Accessories**



17652



17652 removed



**Removable switch mounting**

Adjustable locking stop collar, non-marking type thumb screw lock, adjustable angle two hole switch mount, 3" x 11" aluminum L-rod for adjustment

ITEM #	DESCRIPTION	MOUNT SIZE
17650	REMOVABLE LOCKING SWITCH MOUNT	3/4"
17652	REMOVABLE LOCKING SWITCH MOUNT	7/8"
17654	REMOVABLE LOCKING SWITCH MOUNT	1"




36110



36110  
with sip tube and  
microphone installed  
(not included)



36120



36120  
with sip tube  
installed  
(not included)



36130



36130  
with microphone  
installed  
(not included)

**Microphone and sip/puff mounts**

Durable mounts for ECU's and Augmentative Communication device microphones and sip/puff tubes. Sip/puff mounts are bright nickel plated to prevent rust, have a 3/16" OD tube for vinyl tubing on one end and a 1/8-27 NPT threaded hole on the other end. Mounts to standard 3/8" rod, 5/8" size also available for the sip/puff tube mount. Custom sizes and collar mounts available.

ITEM #	DESCRIPTION	MOUNT COLLAR SIZE	MIC COLLAR SIZE
36110	MICROPHONE AND SIP/PUFF TUBE COMBO MOUNT	3/8"	5/16"
36120	SIP/PUFF TUBE MOUNT	3/8"	
36124	SIP/PUFF TUBE MOUNT	5/8"	
36130	MICROPHONE MOUNT	3/8"	5/16"



**Augmentative Communication Mounting Accessories**



36284




36294

**Double collars**

Fixed double collar mounts for 3/8" rod. The 7/8" size can be combined with mounts in the hardware section (page 41) for various mounting positions on the wheelchair or bed. Also available as a custom double collar with any length spacer between the collars.

ITEM #	DESCRIPTION	COLLAR SIZES
36280	DOUBLE COLLAR - PERPENDICULAR	3/8" - 7/16"
36282	DOUBLE COLLAR - PERPENDICULAR	3/8" - 3/4"
36284	DOUBLE COLLAR - PERPENDICULAR	3/8" - 7/8"
36286	DOUBLE COLLAR - PERPENDICULAR	3/8" - 1"
36290	DOUBLE COLLAR - PARALLEL	3/8" - 7/16"
36292	DOUBLE COLLAR - PARALLEL	3/8" - 3/4"
36294	DOUBLE COLLAR - PARALLEL	3/8" - 7/8"
36296	DOUBLE COLLAR - PARALLEL	3/8" - 1"



34080



Height adjustment for 3/8" rod



34080 mounted on chair


**Removable surface mount combination for 3/8" rod**

Adjustable angle quick release mount for 3/8" rod, 8" long receiving tube allows height adjustment for the 3/8" rod. Surface mount tube is 7/8" diameter and 6" long. Perfect mount for hospital beds and wheelchair rectangular tube seat or back frames.

ITEM #	DESCRIPTION
34080	REMOVABLE SURFACE MOUNT COMBO FOR 3/8" ROD



**Augmentative Communication Mounting Accessories**



34084

Height adjustment for 3/8" rod

34084 with a lollipop mount (see page 40) for wheelchair frame

**Removable tube mount for 3/8" rod**

Adjustable angle quick release mount for 3/8" rod, 8" long receiving tube allows height adjustment for the 3/8" rod. Mounts to wheelchair round 7/8" tube seat or can be combined with mounts in the hardware section (page 39) for various mounting positions on the wheelchair or bed.

ITEM #	DESCRIPTION	MAIN MOUNT COLLAR SIZE
34084	REMOVABLE TUBE MOUNT FOR 3/8" ROD	7/8"



32051



32060



34086



Height adjustment for 3/8" rod

**Removable tube mount for 3/8" rod - separate items**

Create your own mounting or duplicate receivers on a bed frame and a wheelchair, or on a manual and power wheelchair etc.

ITEM #	DESCRIPTION	7/8" TUBE LENGTH
32051	SURFACE MOUNT WITH TUBE - 2 HOLE 3/16" THICK PLATE	2"
32052	SURFACE MOUNT WITH TUBE - 2 HOLE 3/16" THICK PLATE	4"
32050	SURFACE MOUNT WITH TUBE - 2 HOLE 3/16" THICK PLATE	6"
32060	QUICK RELEASE RECEIVER FOR 7/8" TUBE WITH 7/8" COLLAR	
34086	RECEIVING TUBE FOR 3/8" ROD - 8" LONG	



**Augmentative Communication Mounting Accessories**



34088



34088 removed



34088 with a lollipop mount (see page 40) mounted to wheelchair back frame

**Removable small receiver mount for 3/8" rod**

Adjustable angle quick release mount for 3/8" rod, 2" long tube allows height adjustment for the 3/8" rod. Mounts to wheelchair round 7/8" tube seat frames or can be combined with mounts in the hardware section (page 39) for various mounting positions on the wheelchair or bed.

ITEM #	DESCRIPTION	MAIN MOUNT COLLAR SIZE
34088	REMOVABLE SMALL RECEIVER MOUNT FOR 3/8" ROD	7/8"



36314



36314 removed




Adjustable depth for 3/8" rod

**Double collar - removable 3/8"**

Double collar removable mount for 3/8" rod.

ITEM #	DESCRIPTION	MAIN MOUNT COLLAR SIZE
36310	REMOVABLE DOUBLE COLLAR FOR 3/8" ROD	7/16"
36312	REMOVABLE DOUBLE COLLAR FOR 3/8" ROD	3/4"
36314	REMOVABLE DOUBLE COLLAR FOR 3/8" ROD	7/8"
36316	REMOVABLE DOUBLE COLLAR FOR 3/8" ROD	1"



**Augmentative Communication Mounting Accessories**



36252



36252 with gooseneck (not included)




36252 with a lollipop mount (see page 40) mounted to wheelchair back frame

**Gooseneck mount**

Gooseneck mount has a 1/8-27 NPT threaded hole with a 1/4" through hole on the other end. Mounts to wheelchair frame or use a 7/8" size combined with mounts in the hardware section (page 39) for various mounting positions on the wheelchair or bed.

ITEM #	DESCRIPTION	MOUNT SIZE
36250	GOOSENECK MOUNT WITH COLLAR	3/4"
36252	GOOSENECK MOUNT WITH COLLAR	7/8"
36254	GOOSENECK MOUNT WITH COLLAR	1"



36260



36262

**Accessory surface mount with 3/8" collar**

Flat plate mount for batteries, switches, etc. Will accept Velcro for easy attachment or use holes for permanent mounting. Custom plate sizes and hole patterns available.

ITEM #	DESCRIPTION	MOUNTING PLATE SIZE	HOLE CENTER TO CENTER
36260	ACCESSORY SURFACE MOUNT WITH 3/8" COLLAR - SMALL	1" X 2"	1-1/2"
36262	ACCESSORY SURFACE MOUNT WITH 3/8" COLLAR - LARGE	1-1/2" X 3"	2"



36270



36270 mounted on 3/8" rod (not included)

**Gusset support for 3/8" rod**

Double collar support to stabilize 90 degree bends in 3/8" rod, 4" long center rod.

ITEM #	DESCRIPTION
36270	GUSSET SUPPORT FOR 3/8" ROD